

*“Compact Design,
Powerful Results”*

Portable Protection

The Biomist Micro System is a perfect way to sanitize hard-to-reach areas. It can either be pulled by hand or worn backpack-style to allow you room to maneuver. No laborious hand pumping, electricity or batteries are required for operation. Just point and spray! The solution self-dries rapidly so no wiping, rinsing or residual cleanup is necessary.

Thorough, Effortless Coverage

Biomist Sanitizing Systems spray at a constant pressure for precise application, a major advantage over conventional sprayers. Constant pressure generates the atomization of Biomist Formula D2 for even coverage; only the minimally effective amount of sanitizer is used to get the job done. Savings are achieved through reduced time, labor and chemical consumption.

Fast and Effective

Biomist Formula D2 is non-corrosive and safe for use on electronics, metals, and food contact surfaces. Equipment is thoroughly sanitized in seconds without risking damage to sensitive electrical components and wiring. Downtime from water and chemical damage is greatly reduced.

Biomist saves time and eliminates the risk of cross-contamination from wiping. Equipment and portable items that carry pathogens from one point to another can be sanitized quickly and efficiently.



Saving Money, Saving Lives.

Biomist[®] Micro Sanitizing System



Protect your Brand and your Business

Emerging regulations place greater accountability on companies to protect the public, as uncontrolled pathogens cost billions of dollars and tens of thousands of lives per year. Current sanitizing practices are costly, time-consuming, often ineffective, and sometimes dangerous. Biomist strengthens sanitation programs and helps you maintain a cleaner and safer environment. Use Biomist proactively to guard your reputation and minimize the risk associated with recalls and negative media.

Affordable Assurance

It is well known that bacteria and viruses can survive on surfaces for extended periods of time and Biomist is the fastest sanitizing method available. It provides the best of both worlds: Biomist is more effective and less expensive to use than all previous sanitizing methods while being safe for employees, equipment and the environment.

Tested and Approved

The Biomist Micro Sanitizing System is registered with the FDA as a medical device. It is used by hospitals, pharmaceutical plants, and several of the nation's largest food processors. Independent testing has proved it effective in killing:

Staphylococcus aureus MRSA (ATCC 33591)
Staphylococcus aureus (ATCC 6538)
Enterococcus faecium VRE (ATCC 51559)
Escherichia coli 0157:H7 (ATCC 43895)
Pseudomonas aeruginosa (ATCC 15442)
Listeria monocytogenes (ATCC 19111)
Salmonella choleraesuis (ATCC 10708)
HIV- 1 (associated with AIDS)
Hepatitis B Virus (HBV)
Norovirus
and the vast majority of all other viruses and bacteria

The Biomist Micro Sanitizing System is certified compliant with most relevant industry standards.

Kill Germs Fast

Biomist Formula D2 contains a concentrated solution of alcohol (long recognized as an optimal sanitizing agent), as well as a four-chain quaternary ammonium compound to continue sanitizing after the alcohol has evaporated.

Fire-Safe Operation

To eliminate the problem of alcohol's flammability, Biomist's patented technology uses high pressure bottled carbon dioxide (CO2) as a carrier to spray the sanitizer. By enveloping the alcohol mist in a stream of CO2 propellant, Biomist displaces the oxygen needed for combustion. It can be used in the presence of open flame without danger.

Use It Everywhere

- Health Care Institutions
- Travel and Hospitality
- Food Plants
- Ambulance Services
- Clean Rooms and Laboratories
- Restaurants
- and other critical environments



The fastest, most effective sanitizing method available.



Biomist, Inc.
573 North Wolf Road
Wheeling, IL 60090

Phone: 847-850-5530

Fax: 847-850-5535

E-mail: info@biomistinc.com

Or log on to www.biomistinc.com

BIOMIST Micro

Saving Money, Saving Lives.